réf. CCC\_35

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## Directions for use Prosthetic torque wrench

## 1 - DESCRIPTION

Made of stainless steel which can be sterilized. Resistance enables an optimal protection against the aggression of the decontamination and of the sterilization products.

Dismantled, it can be easily cleaned

Several torque adjustments are available : 10N.cm, 15N.cm, 20N.cm, 25N.cm, 30N.cm, 35N.cm, 40N.cm.

Tolerance applied to the torque value is more or less 2 N.cm.

Instrument approved for 5000 uses

## 2 - INSTRUCTIONS FOR USE

Calibrated tightening of the prosthetic parts. If the assembly of prosthetic constructions are done without the recommended prestressing (cf instructions for use), the screw could get unscrewed after a while. Moreover, a too weak torque would damage the screw pitch and would not maintain the prosthetic part, a too heavy prevailing torque would not allow to absorb the chewing of stress.

The torque adjustment has to be done by turning the knurling wheel **7** to show the desired torque value and to make the straight line corresponding to this value match printed of the torque wrench body.

The mark arrow located on the head of the wrench **1** gives the screwing direction.

The torque wrench has to be calibrated once a year to be precisely controlled. To make it calibrated, please send it to the head office of the company.

The torque wrench is guaranteed for one year for parts and labour against any manufacturing defect (not for use defect) from the delivery date.

Storage : When you don't use the torque wrench, adjust it on the minimum torque value to avoid any modification of the elasticity properties of the spring and which could cause any torque variation.

The key bends when the desired tightening torque is reached. You should not apply stress any more on the handle at this time.

## 3 - MAINTENANCE

The torque wrench is sold non-sterile

It has to be cleaned up, disinfected and sterilized before use in accordance with the directions for use below.

### 3-1 Cleaning

Immediately after using the wrench, put it into a container filled up with a disinfectant/de-

tergent solution and put a linen on it. The aim is to prevent from the drying of the infectious element from the patient, to dissolve it and make the cleaning easier. The use of chlorine-based cleaning or disinfecting products (bleach) has to be avoided.

### 3-2 Disinfection: Choose between two solutions (suggested solutions)

To carry out a thermic disinfection, contact with boiling water (more than 90°) during 3 to 10 minutes. Then, rinse all the parts with cold distilled water.

To carry out a chemical disinfection with an appropriate product, which does not damage steel (fungicide and virucidal) acting quickly. Then rinse all the parts with cold distilled water.

Completely dismantle the wrench taking care to respect the dismantle order described in chapter 4.

Clean up all the external and internal surfaces of the instrument with a soft bristle brush under hot water.

Rinse all the less accessible holes of the head **1** by injecting some hot water with a syringe without needle. If necessary, do the same for the inner parts of the shaft **5** and of the knurling wheel **7**.

If the notch has been dismantled, clean up the wrench by ultrasound with a disinfectant/detergent solution appropriate for stainless steel. While the parts are cleaned up, avoid the contact between each of them.

Rinse all the parts with demineralized water to eliminate all the detergents and particles traces.

Remark : All blood and layers residues make the disinfection and sterilization less efficient. Avoid any contact between the instrument and other nickelled instruments. During all the cleaning cycle avoid any projection of liquid and work with the appropriate protections.

Be Careful : If the torque wrench is not perfectly cleaned, particularly in the area situated on the shank between the body and the head of the torque wrench, we can't guarantee the values given by the wrench. Without a regular cleaning, the value can increase.

### 3-2 Preparation before sterilization

Important: All the parts have to be carefully and meticulously dried before sterilization, never reassemble a torque wrench while its components are not perfectly dry.

Carefully reassemble the torque wrench carefully respecting the indication of chapter 5.

Adjust the wrench on the minimum torque value.

Carefully inspect the wrench to be sure there is no dirty parts left.

Pack the wrench.

### 3-3 Sterilization

The etK torque wrench has to be sterilized in autoclave before use (during 18 minutes at 134°C (273,2°F)). Sterilization in accordance with the sterilization cycles given by the autoclave manufacturer. The sterilization process must have been validated.

We recommend to use devices with evacuating pump in order to reduce the risk of air areas formation. It is particularly important for all the hollow instruments and to guarantee a perfect dry. If some corrosion traces appear, all the components have to be treated with a 0,1% sodium nitrate solution before the sterilization process.

The sterilization with autoclave is the only method recommended by the French Ministry of Health circulars DGS/5C/DHOS/E2/2001 N°138 of March, 14th, 2001.

Remark: The sterilizer with hot air is forbidden because it accelerates the ageing process of the spring and can provoke a modification of the real torque applied.

## 4 - DISMANTLING OF THE TORQUE WRENCH

**Step 1** : Unscrew the knurling wheel **7** and remove the adjustment spring **6** from the torque wrench.

**Step 2** : Unscrew the assembling axle **2** (with an external hexagonal screwdriver **9** type CCL\_HE12.xx).

**Step 3** : Remove the Rod **5** and the adjusment spring **4**.

Be Careful : Dismantling is necessary only if the pawl is jammed.

The reassembly is delicate, avoid to dismantle this part if the notch normally function.

## 5 - ASSEMBLY OF THE TORQUE WRENCH

Be Careful : It is necessary to use oil for contra-angle and pieces on the indicated points by **8**.

**Step 1** : Assemble the whole spring **4** and notch **3** in the head of the wrench, (**be careful, the notch can only be assemble on one way**. In order to facilitate the assembly, use a hexagonal screwdriver.

**Step 2** : Put the shaft **1** in the head of the wrench **1** and screw the assembling axle **2** Screw without forcing. (20 N.cm)

Be Careful : Hand screwing with a external hexagonal screwdriver **9**, never use a pawl wrench to screw this axle.

## FUNCTIONING CONTROL

A checking test of the assembling and of the functioning is necessary before the surgery, before sterilization, and while the instruments are prepared.

When the torque wrench has been assembled and adjusted on several torque values between 10 and 40N.cm, check by, while holding the head of the wrench in your hand and applying a perpendicular stress to the axle of the wrench on the knurling wheel that the torque wrench sets off and is not blocked.